

# Anti-RAIDD Rabbit antibody



Catalog Number: 167815

## Product name

Anti-RAIDD Rabbit antibody

## Specificity

Human

## Antibody description

RAIDD Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human RAIDD

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

IgG

## Applications

WB, IHC-P, ICC/IF, FC, IP

## Dilutions

WB: 1/1000


IHC: 1/20

ICC/IF: 1/50

FC: 1/20

IP: 1/20

## Validations



Western blot detection of RAIDD in HeLa, A549, MCF-7, 293T cell lysates using RAIDD Rabbit mAb (1:1000 diluted). Predicted band size: 23 kDa. Observed band size: 23 kDa.