

# Anti-MLLT6 Rabbit antibody



Catalog Number: 167817

## Product name

Anti-MLLT6 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

MLLT6 Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human MLLT6

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of MLLT6 in HeLa,A549,HL-60,U2OS cell lysates using MLLT6 Rabbit mAb(1:1000 diluted).Predicted band size:112kDa.Observed band size:112kDa.