

Anti-FRS2 Rabbit antibody



Catalog Number: 167824

Product name

Anti-FRS2 Rabbit antibody

Specificity

Human

Antibody description

FRS2 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human FRS2

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, FC, IP


Dilutions

WB: 1/1000

FC: 1/20

IP: 1/20

Validations



Western blot detection of FRS2 in HeLa cell lysates using FRS2 Rabbit mAb(1:1000 diluted). Predicted band size:57kDa.Observed band size:57kDa.