

Anti-UBA3 Rabbit antibody



Catalog Number: 167832

Product name

Anti-UBA3 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

UBA3 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human UBA3

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/2000-1/10000

Validations



Western blot detection of UBA3 in HeLa,A549,HL-60 cell lysates using UBA3 Rabbit mAb(1:1000 diluted).Predicted band size:52kDa.Observed band size:52kDa.