

# Anti-Nicotinic Acetylcholine Receptor alpha 4/CHRNA4 Rabbit antibody



Catalog Number: 167833

## Product name

Rabbit IgG

Anti-Nicotinic Acetylcholine Receptor alpha 4/CHRNA4 Rabbit antibody

## Applications

WB, ICC/IF, IP

## Specificity

Human, Mouse, Rat

## Dilutions

WB: 1/2000-1/10000

## Antibody description

Nicotinic Acetylcholine Receptor alpha 4/CHRNA4 Rabbit monoclonal antibody


ICC/IF: 1/20

IP: 1/20

## Preparation

Antigen: A synthetic peptide of human Nicotinic Acetylcholine Receptor alpha 4/CHRNA4

## Validations



## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Western blot detection of Nicotinic Acetylcholine Receptor alpha 4/CHRNA4 in HeLa, CHO-K1, C6, Ramos cell lysates using Nicotinic Acetylcholine Receptor alpha 4/CHRNA4 Rabbit mAb(1:1000 diluted). Predicted band size:70kDa. Observed band size:70kDa.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type