

Anti-phospho-Rb (Thr356) Rabbit antibody



Catalog Number: 167841

Product name

Anti-phospho-Rb (Thr356) Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human

Applications

WB, IHC-P, IP

Antibody description

phospho-Rb (Thr356) Rabbit monoclonal antibody

Dilutions

WB: 1/1000

IHC: 1/50

IP: 1/20

Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Thr356 of human Rb

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Validations



Western blot detection of Rb (phospho-Thr356) in HL-60, U2OS, C6 using Rb (phospho-Thr356) Rabbit mAb(1:1000 diluted)