

Anti-Annexin-4/ANXA4 Rabbit antibody



Catalog Number: 167846

Product name

Rabbit IgG

Anti-Annexin-4/ANXA4 Rabbit antibody

Applications

Specificity

WB, IHC-P, ICC/IF, FC

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/1000

Annexin-4/ANXA4 Rabbit monoclonal antibody

IHC: 1/500

Preparation

ICC/IF: 1/20

Antigen: Recombinant protein of human Annexin-4/ANXA4

FC: 1/50

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage


Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Validations



Western blot detection of Annexin A4 in HeLa,CHO-K1,C6 using Annexin A4 Rabbit mAb(1:1000 diluted)