

Anti-MST1 Rabbit antibody



Catalog Number: 167847

Product name

Anti-MST1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

MST1 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human MST1

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, FC


Dilutions

WB: 1/1000-1/5000

IHC: 1/20-1/50

FC: 1/20

Validations



Western blot detection of MST1 in HeLa,C6,Rat Brain lysates using MST1 Rabbit mAb(1:1000 diluted).Predicted band size:80kDa.Observed band size:80kDa.