

Anti-Rad50 Rabbit antibody



Catalog Number: 167858

Product name

Anti-Rad50 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Rad50 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human Rad50

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

IgG

Applications


WB, IHC-P

Dilutions

WB: 1/1000

IHC: 1/500

Validations



Western blot detection of Rad50 in HeLa,A549,HL-60 using Rad50 Rabbit mAb(1:1000 diluted)