

# Anti-Glucose 6 Phosphate Dehydrogenase Rabbit antibody



Catalog Number: 167874

## Product name

IgG

Anti-Glucose 6 Phosphate Dehydrogenase Rabbit antibody

## Applications

WB, IHC-P, ICC/IF, FC

## Specificity

Human

## Dilutions

WB: 1/1000

## Antibody description

Glucose 6 Phosphate Dehydrogenase Rabbit monoclonal antibody

IHC: 1/20

ICC/IF: 1/20

## Preparation

Antigen: A synthetic peptide of human Glucose 6 Phosphate Dehydrogenase

FC: 1/20

## Validations

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA


## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type



Western blot detection of Glucose 6 Phosphate Dehydrogenase in HeLa, A549, U2OS using Glucose 6 Phosphate Dehydrogenase Rabbit mAb(1:1000 diluted)