

# Anti-Arp2 Rabbit antibody



Catalog Number: 167883

## Product name

Anti-Arp2 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Arp2 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human Arp2

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

IgG

## Applications

WB, IHC-P, ICC/IF


## Dilutions

WB: 1/1000

IHC: 1/50

ICC/IF: 1/50

## Validations



Western blot detection of Arp2 in C6, mouse Brain, Ramos using Arp2 Rabbit mAb(1:1000 diluted)