

# Anti-Eph receptor B4/HTK Rabbit antibody



Catalog Number: 167887

## Product name

Anti-Eph receptor B4/HTK Rabbit antibody

## Specificity

Human, Mouse

## Antibody description

Eph receptor B4/HTK Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human Eph receptor B4/HTK

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC-P

## Dilutions

WB: 1/1000

IHC: 1/100

## Validations



Western blot detection of Eph receptor B4/HTK in HeLa,A549,C6 cell lysates using Eph receptor B4/HTK Rabbit mAb(1:1000 diluted).Predicted band size:108kDa.Observed band size:135kDa.