

# Anti-LYRIC/AEG1 Rabbit antibody



Catalog Number: 167897

## Product name

Anti-LYRIC/AEG1 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

LYRIC/AEG1 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human LYRIC/AEG1

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

IgG

## Applications

WB, IHC-P, ICC/IF


## Dilutions

WB: 1/2000-1/10000

IHC: 1/50

ICC/IF: 1/100

## Validations



Western blot detection of LYRIC/AEG1 in CHO-K1, Ramos using LYRIC/AEG1 Rabbit mAb(1:1000 diluted)