

# Anti-TFPI Rabbit antibody



Catalog Number: 167907

## Product name

Anti-TFPI Rabbit antibody

## Specificity

Human

## Antibody description

TFPI Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human TFPI

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of TFPI in A549, MCF-7 cell lysates using TFPI Rabbit mAb(1:1000 diluted). Predicted band size: 35kDa. Observed band size: 40-55kDa.