

Anti-RNA polymerase II CTD repeat YSPTSPS Rabbit antibody



Catalog Number: 167914

Product name

Anti-RNA polymerase II CTD repeat YSPTSPS Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

RNA polymerase II CTD repeat YSPTSPS Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human RNA polymerase II CTD repeat YSPTSPS

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/1000


IHC: 1/100

ICC/IF: 1/100

FC: 1/20

IP: 1/20

Validations



Western blot detection of RNA polymerase II CTD repeat YSPTSPS in HeLa, A549, MCF-7, U251, HL-60 cell lysates using RNA polymerase II CTD repeat YSPTSPS Rabbit mAb (1:1000 diluted). Predicted band size: 192 kDa. Observed band size: 250 kDa.