

Anti-HLF Rabbit antibody



Catalog Number: 167938

Product name

Anti-HLF Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

HLF Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human HLF

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of HLF in HeLa,A549,HL-60,U87-MG cell lysates using HLF Rabbit mAb(1:1000 diluted).Predicted band size:33kDa.Observed band size:33kDa.