

Anti-Maltose Binding Protein Rabbit antibody



Catalog Number: 167939

Product name

Monoclonal

Anti-Maltose Binding Protein Rabbit antibody

Ig Type

Specificity

Rabbit IgG

Escherichia Coli

Applications

Antibody description

WB

Maltose Binding Protein Rabbit monoclonal antibody


Dilutions

WB: 1/500-1/2000

Preparation

Validations

Antigen: A synthesized peptide derived from human Maltose Binding Protein



Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of Maltose Binding Protein in HeLa,A549 cell lysates using Maltose Binding Protein Rabbit mAb(1:1000 diluted).Observed band size:43kDa.

Clonality