

# Anti-TORC2 Rabbit antibody



Catalog Number: 167940

## Product name

Anti-TORC2 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

TORC2 Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human TORC2

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of TORC2 in HeLa,A549,HL-60,U251 cell lysates using TORC2 Rabbit mAb(1:1000 diluted).Predicted band size:73kDa.Observed band size:80kDa.