

# Anti-beta Dystroglycan Rabbit antibody



Catalog Number: 167944

## Product name

Anti-beta Dystroglycan Rabbit antibody

## Specificity

Human

## Antibody description

beta Dystroglycan Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human beta Dystroglycan

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of beta Dystroglycan in MCF7, HepG2, Lncap cell lysates using beta Dystroglycan Rabbit mAb(1:1000 diluted). Predicted band size: 27kDa. Observed band size: 42-44kDa.