

# Anti-Aconitase 2 Rabbit antibody



Catalog Number: 167972

## Product name

Anti-Aconitase 2 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Aconitase 2 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human Aconitase 2

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

IgG

## Applications

WB, IHC-P, ICC/IF, FC

## Dilutions


WB: 1/2000-1/10000

IHC: 1/50

ICC/IF: 1/20

FC: 1/20

## Validations



Western blot detection of Aconitase 2 in HL-60,C6 using Aconitase 2 Rabbit mAb(1:1000 diluted)