

Anti-NDUFS2 Rabbit antibody



Catalog Number: 167975

Product name

Anti-NDUFS2 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

NDUFS2 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human NDUFS2

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IP


Dilutions

WB: 1/1000

IHC: 1/20-1/100

IP: 1/20

Validations



Western blot detection of NDUFS2 in HeLa, HepG2, U2OS, U251, HL-60 using NDUFS2 Rabbit mAb(1:1000 diluted)