

# Anti-Ube2L3/UBCH7 Rabbit antibody



Catalog Number: 167978

## Product name

Rabbit IgG

Anti-Ube2L3/UBCH7 Rabbit antibody

## Applications

## Specificity

WB, IHC-P, ICC/IF, IP

Human, Mouse, Rat

## Dilutions

## Antibody description

WB: 1/1000

Ube2L3/UBCH7 Rabbit monoclonal antibody

IHC: 1/50

## Preparation

ICC/IF: 1/50

Antigen: A synthetic peptide of human Ube2L3/UBCH7

IP: 1/20

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage


Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

## Validations



Western blot detection of Ube2L3/UBCH7 in HeLa, A549, HL-60, U2OS, C6 using Ube2L3/UBCH7 Rabbit mAb(1:1000 diluted)