

Anti-PNK/PNKP Rabbit antibody



Catalog Number: 168002

Product name

Anti-PNK/PNKP Rabbit antibody

Specificity

Human, Mouse

Antibody description

PNK/PNKP Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human PNK/PNKP

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IP


Dilutions

WB: 1/1000

IHC: 1/20

IP: 1/20

Validations



Western blot detection of PNK/PNKP in HeLa,A549,HL-60,U251 cell lysates using PNK/PNKP Rabbit mAb(1:1000 diluted).Predicted band size:57kDa.Observed band size:57kDa.