

Anti-LTK Rabbit antibody



Catalog Number: 168007

Product name

Anti-LTK Rabbit antibody

Specificity

Human

Antibody description

LTK Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human LTK

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

WB: 1/1000

IHC: 1/50

Validations



Western blot detection of LTK in A549,HL-60,U2OS,C6 cell lysates using LTK Rabbit mAb(1:1000 diluted).Predicted band size:92kDa.Observed band size:92kDa.