

Anti-PBK/SPK Rabbit antibody



Catalog Number: 168011

Product name

Anti-PBK/SPK Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

PBK/SPK Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human PBK/SPK

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

IgG

Applications

WB, IHC-P, FC, IP

Dilutions


WB: 1/1000

IHC: 1/500

FC: 1/50

IP: 1/20

Validations



Western blot detection of PBK/SPK in HeLa,A549,HL-60,U2OS cell lysates using PBK/SPK Rabbit mAb(1:1000 diluted).Predicted band size:36kDa.Observed band size:40kDa.