

# Anti-HIF Prolyl Hydroxylases Rabbit antibody



Catalog Number: 168024

## Product name

Rabbit IgG

Anti-HIF Prolyl Hydroxylases Rabbit antibody

## Applications

## Specificity

WB, IHC-P, FC, IP

Human

## Dilutions

## Antibody description

WB: 1/1000

HIF Prolyl Hydroxylases Rabbit monoclonal antibody

IHC: 1/50

FC: 1/50

## Preparation

IP: 1/20

Antigen: A synthetic peptide of human HIF Prolyl Hydroxylases

## Validations

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA


## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type



Western blot detection of HIF Prolyl Hydroxylases in HeLa,A549,HL-60,U251,U87-MG cell lysates using HIF Prolyl Hydroxylases Rabbit mAb(1:1000 diluted).Predicted band size:57kDa.Observed band size:43kDa.