

# Anti-Lgi1/EPT Rabbit antibody



Catalog Number: 168026

## Product name

Rabbit IgG

Anti-Lgi1/EPT Rabbit antibody

## Applications

## Specificity

WB, IHC-P, IP

Human, Mouse, Rat

## Dilutions

## Antibody description

WB: 1/1000

Lgi1/EPT Rabbit monoclonal antibody

IHC: 1/50

## Preparation

IP: 1/20

Antigen: A synthetic peptide of human Lgi1/EPT

## Validations

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA


## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type



Western blot detection of Lgi1/EPT in HeLa,A549 cell lysates using Lgi1/EPT Rabbit mAb(1:1000 diluted).Predicted band size:64kDa.Observed band size:64kDa.