

Anti-PHD1/prolyl hydroxylase Rabbit antibody



Catalog Number: 168031

Product name

Rabbit IgG

Anti-PHD1/prolyl hydroxylase Rabbit antibody

Applications

Specificity

WB, IHC-P, ICC/IF, FC

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/500-1/2000

PHD1/prolyl hydroxylase Rabbit monoclonal antibody

IHC: 1/50-1/200

Preparation

ICC/IF: 1/50-1/200

Antigen: A synthesized peptide derived from human PHD1/prolyl hydroxylase

FC: 1/50

Validations

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of PHD1/prolyl hydroxylase in HeLa, A549, HL-60, U2OS, U87-MG cell lysates using PHD1/prolyl hydroxylase Rabbit mAb (1:1000 diluted). Observed band size: 44 kDa.