

# Anti-PHD1/prolyl hydroxylase Rabbit antibody



Catalog Number: 168031

## Product name

IgG

Anti-PHD1/prolyl hydroxylase Rabbit antibody

## Applications

## Specificity

WB, IHC-P, ICC/IF, FC

Human, Mouse, Rat

## Dilutions

## Antibody description

WB: 1/500-1/2000

PHD1/prolyl hydroxylase Rabbit monoclonal antibody

IHC: 1/50-1/200

ICC/IF: 1/50-1/200

## Preparation

FC: 1/50

Antigen: A synthesized peptide derived from human PHD1/prolyl hydroxylase

## Validations

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.


## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type



Western blot detection of PHD1/prolyl hydroxylase in HeLa, A549, HL-60, U2OS, U87-MG cell lysates using PHD1/prolyl hydroxylase Rabbit mAb (1:1000 diluted). Observed band size: 44 kDa.