

Anti-Cytokeratin 9 Rabbit antibody



Catalog Number: 168046

Product name

Anti-Cytokeratin 9 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Cytokeratin 9 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human Cytokeratin 9

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IP


Dilutions

WB: 1/1000

IHC: 1/50

IP: 1/20

Validations



Western blot detection of Cytokeratin 9 in C6,Ramos using Cytokeratin 9 Rabbit mAb(1:1000 diluted)