

Anti-DAB2 Rabbit antibody



Catalog Number: 168047

Product name

Anti-DAB2 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

DAB2 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human DAB2

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

IgG

Applications

WB, IHC-P, IP


Dilutions

WB: 1/1000

IHC: 1/50

IP: 1/20

Validations



Western blot detection of DAB2 in
HepG2,U2OS,U251 using DAB2 Rabbit
mAb(1:1000 diluted)