

Anti-AK2 Rabbit antibody



Catalog Number: 168049

Product name

Anti-AK2 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

AK2 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human AK2

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, IP

Dilutions


WB: 1/1000

IHC: 1/50

ICC/IF: 1/100

IP: 1/20

Validations



Western blot detection of AK2 in HeLa,CHO-K1,C6,mouse Brain,Ramos using AK2 Rabbit mAb(1:1000 diluted)