

Anti-CKMT1 (1A6) Mouse antibody



Catalog Number: 168069

Product name

Anti-CKMT1 (1A6) Mouse antibody

Specificity

Human, Mouse

Antibody description

CKMT1 (1A6) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human CKMT1 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of CKMT1 in Mouse Colon, 293 and A431 cell lysates using CKMT1 mouse mAb (1:1000 diluted). Predicted band size: 47KDa. Observed band size: 47KDa.