

# Anti-Hexokinase II (1E8) Mouse antibody



Catalog Number: 168070

## Product name

Anti-Hexokinase II (1E8) Mouse antibody

## Specificity

Rat, Mouse, Monkey

## Antibody description

Hexokinase II (1E8) Mouse monoclonal antibody

## Preparation

Antigen: Recombinant human Hexokinase II protein.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

WB

## Dilutions

WB: 1/500

## Validations



Western blot detection of Hexokinase II in 3T3, COS7 and C6 cell lysates using Hexokinase II mouse mAb(dilution 1:500). Predicted band size: 102kDa. Observed band size: 102kDa.