

Anti-C Reactive Protein (1H3) Mouse antibody



Catalog Number: 168074

Product name

Anti-C Reactive Protein (1H3) Mouse antibody

Specificity

Human

Antibody description

C Reactive Protein (1H3) Mouse monoclonal antibody

Preparation

Antigen: Native C-Reactive Protein.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300.

Storage

Store at 4°C short term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG1


Applications

ELISA

Validations



sELISA standard Curve for C-Reactive Protein: Capture Antibody Mouse mAb to C-Reactive Protein at 1ug/ml and 168074 was used for detecting.



sELISA standard Curve for C-Reactive Protein: Capture Antibody Rabbit pAb to C-Reactive Protein at 4ug/ml and 168074 was used for detecting.