

Anti-BRAP (1E7) Mouse antibody



Catalog Number: 168083

Product name

Anti-BRAP (1E7) Mouse antibody

Specificity

Human

Antibody description

BRAP (1E7) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human BRAP protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

IgG1


Applications

WB


Dilutions

WB: 1/1000;Other 1/500

Validations



Western blot detection of BRAP in A549 and HUVEC cell lysates using BRAP mouse mAb (1:500 diluted). Predicted band size: 67KDa. Observed band size: 67KDa.



Western blot analysis of extracts from CHO-K1 cells, transfected with a human pFLAG-CMV2-BRAP construct (A) or transfected with a human pFLAG-CMV2 construct (B), using BRAP mouse mAb (1:1000 diluted).