

# Anti-JMJD6 (3G5) Mouse antibody



Catalog Number: 168090

## Product name

Anti-JMJD6 (3G5) Mouse antibody

## Specificity

Human

## Antibody description

JMJD6 (3G5) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human JMJD6(N-terminus) fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


WB, ICC/IF

## Dilutions


WB: 1/1000

;ICC: 1/200

## Validations



Immunocytochemistry of HeLa cells using anti-JMJD6(N-terminus) mouse mAb diluted 1:200.



Western blot detection of JMJD6(N-terminus) in HeLa and HL-60 lysates using JMJD6(N-terminus) mouse mAb (1:1000 diluted). Predicted band size: 46KDa. Observed band size: 62KDa.