

Anti-SMC1A (4C5) Mouse antibody



Catalog Number: 168092

Product name

Anti-SMC1A (4C5) Mouse antibody

Specificity

Human

Antibody description

SMC1A (4C5) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human SMC1A(C-term.) protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


WB, ICC/IF

Dilutions


WB: 1/1000

;ICC: 1/100

Validations



Western blot detection of SMC1A(C-term) in MOLT-4,U87 MG and HeLaNE cell lysates using SMC1A (N-terminus) mouse mAb (1:1000 diluted).Predicted band size: 143KDa.Observed band size: 143KDa.



Immunocytochemistry staining of HeLa cells fixed with 4% Paraformaldehyde and using anti-SMC1A mouse mAb(dilution 1:100).