

Anti-CBR1 (2C9) Mouse antibody



Catalog Number: 168095

Product name

Anti-CBR1 (2C9) Mouse antibody

Specificity

Human

Antibody description

CBR1 (2C9) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human CBR1 protein fragments expressed in E.coli

Formulation

Purified mouse monoclonal in PBS(pH 7.4)containing with 0.02% sodium azide 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


WB, ICC/IF

Dilutions


WB: 1/1000

;ICC: 1/100

Validations



Western blot detection of CBR1 in HeLa,A431 and MDA-MB-468 cell lysates using CBR1 mouse mAb (1:1000 diluted).Predicted band size:30KDa,Observed band size:30KDa.



Immunocytochemistry stain of HeLa using CBR1 mouse mAb (1:100).