

# Anti-ABCB5 (8D2) Mouse antibody



Catalog Number: 168098

## Product name

Mouse IgG1

Anti-ABCB5 (8D2) Mouse antibody

## Applications

## Specificity

WB, IHC-P, ICC/IF

Human

## Dilutions

## Antibody description

WB: 1/2000

ABCB5 (8D2) Mouse monoclonal antibody

IF: 1/200

## Preparation

IHC: 1/50-300

Antigen: Synthetic Peptide of ABCB5

## Validations

## Formulation

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.


## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type



Western blot detection of ABCB5 in K562 cell lysates using ABCB5 mouse mAb (1:2000 diluted). Predicted band size: 90kDa. Observed band size: 90kDa.