

# Anti-eIF4A1 (6D7) Mouse antibody



Catalog Number: 168099

## Product name

Anti-eIF4A1 (6D7) Mouse antibody

## Specificity

Human, Mouse, Rat

## Antibody description

eIF4A1 (6D7) Mouse monoclonal antibody

## Preparation

Antigen: Synthetic Peptide of eIF4A1

## Formulation

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

WB, IHC-P, ICC/IF


## Dilutions

WB: 1/1000-3000


IF: 1/100-200

IHC: 1/50-300

## Validations



IF analysis of HeLa with (Left) and DAPI (Right)



diluted at 1:100

Western blot detection of eIF4A1 in T47D, 3T3, COS7, C6 and HeLa cell lysates using eIF4A1 mouse mAb (1:1000 diluted). Predicted band size: 48KDa. Observed band size: 48KDa.