

Anti-Leptin (7E3) Mouse antibody



Catalog Number: 168104

Product name

Anti-Leptin (7E3) Mouse antibody

Specificity

Human

Antibody description

Leptin (7E3) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant full length of human Leptin protein expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

IgG1

Applications

ELISA

Validations



Indirect ELISA assay for anti-Leptin mouse mAb. Antigen coating concentration: 4ug/ml.