

Anti-kan (4E9) Mouse antibody

Catalog Number: 168106

Product name

Anti-kan (4E9) Mouse antibody

Antibody description

kan (4E9) Mouse monoclonal antibody

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of nptII in non-trangenic Brassica napus L. and transgenic Brassica napus L. cell lysates using nptII mouse mAb (1:1000 diluted). Predicted band size: 29KDa. Observed band size: 29KDa.