

Anti-EGFR (1A11) Mouse antibody



Catalog Number: 168108

Product name

Anti-EGFR (1A11) Mouse antibody

Specificity

Human

Antibody description

EGFR (1A11) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human EGFR protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide 0.1%BSA and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB, IP


Dilutions

WB: 1/2000

Validations



Immunoprecipitation analysis of HeLa cell lysates



using EGFR mouse mAb.

Western blot detection of EGFR in A549 cell lysates using EGFR mouse mAb(dilution 1:2000).Predicted band size:134 Kda.Observed band size:175KDa.