

Anti-Human IgG (2D7) Mouse antibody



Catalog Number: 168109

Product name

Anti-Human IgG (2D7) Mouse antibody

Specificity

Human

Antibody description

Human IgG (2D7) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human Human IgG protein fragments expressed in E.coli.

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300.

Storage

Store at 4°C short term.Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG1


Applications

ELISA

Validations



Recombinant Anthrax Protective Antigen used as Coating antigen, followed by different dilutions of supernatants which includes secreted Anti-Anthrax Protective Antigen Human antibody, HRP conjugated Anti-Human IgG(168109) Mouse mAb was used for coloring.



Recombinant Anthrax Protective Antigen used as Coating antigen, followed by Anti-Anthrax Protective Antigen Human antibody and then incubated with Anti-Human IgG(168109) Mouse mAb, HRP conjugated goat anti Mouse IgG was used for coloring.