

# Anti-PI 3 Kinase p85 beta (8D9) Mouse antibody



Catalog Number: 168120

## Product name

Anti-PI 3 Kinase p85 beta (8D9) Mouse antibody

## Ig Type

Mouse IgG1

## Specificity

Human, Mouse, Rat

## Applications

WB

## Antibody description

PI 3 Kinase p85 beta (8D9) Mouse monoclonal antibody

## Dilutions

WB: 1/1000

## Preparation

Antigen: Purified recombinant human PI3 Kinase p85 beta protein fragments expressed in E.coli.

## Validations

## Formulation


PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal



Western blot detection of PI3 Kinase p85 beta in C6,K562,MCF7,U87MG,HCT116,3T3 and HeLa cell lysates using PI3 Kinase p85 beta mouse mAb (1:1000 diluted).Predicted band size:85KDa.Observed band size:85KDa.