

# Anti-IMP3 (1F1) Mouse antibody



Catalog Number: 168121

## Product name

Anti-IMP3 (1F1) Mouse antibody

## Ig Type

IgG1

## Specificity

Human

## Applications

WB

## Antibody description

IMP3 (1F1) Mouse monoclonal antibody

## Dilutions

WB: 1/1000

## Preparation

Antigen: Purified recombinant human IMP3 protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.


## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Validations



Western blot detection of IMP3 in MCF7,SKVO3,NTERA2,Jurkat,Hela and K562 cell lysates using IMP3 mouse mAb (1:1000 diluted).Predicted band size:70KDa.Observed band size:70KDa.