

Anti-HSP90 alpha (6H7) Mouse antibody



Catalog Number: 168123

Product name

Anti-HSP90 alpha (6H7) Mouse antibody

Specificity

Human, Mouse, Rat, Monkey

Antibody description

HSP90 alpha (6H7) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human Hsp90 alpha protein fragments expressed in E.coli

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide 1mg/ml BSA and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

IgG1


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of Hsp90 alpha in R Brain, Jurkat, CHO-K1, COS7, C6, 3T3 and HeLa cell lysates using Hsp90 alpha mouse mAb (1:1000 diluted). Predicted band size: 90KDa. Observed band size: 90KDa.