

Anti-Cytokeratin 16 (10C7) Mouse antibody



Catalog Number: 168134

Product name

Anti-Cytokeratin 16 (10C7) Mouse antibody

Specificity

Human, Mouse, Rat

Antibody description

Cytokeratin 16 (10C7) Mouse monoclonal antibody

Preparation

Antigen: Synthetic Peptide of CK16

Formulation

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


IHC-P, WB

Dilutions


IHC: 1/200

IF: 1/50-200

Validations



IHC staining of human gullet cancer tissue with CK16 mouse mAb(6F6) diluted at 1:200.



Western blot detection of Cytokeratin 16 in T47D and HeLa cell lysates using Cytokeratin 16 mouse mAb (1:3000 diluted). Predicted band size: 51kDa. Observed band size: 51kDa.