

Anti-HDAC1 (3F5) Mouse antibody



Catalog Number: 168147

Product name

Anti-HDAC1 (3F5) Mouse antibody

Specificity

Human, Rat, Mouse

Antibody description

HDAC1 (3F5) Mouse monoclonal antibody

Preparation

Antigen: Synthetic peptide conjugated to KLH.

Formulation

PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

IHC-P


Dilutions

IHC: 1/100-200

Validations



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using HDAC1 Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using HDAC1 Mouse mAb diluted at 1:200.